Otto Bock®

QUALITY FOR LIFE

CE

B 600

 (D)	Bedienungsanleitung	1
 (GB USA)	Instructions for Use	55

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Table of Contents

	Page
1 General Information	57
2 Declaration of Conformity	57
3 Safety Instructions	58
3.1 Warning Symbols	58
3.2 General Safety Instructions.....	58
4 Field of Application	61
5 Service	62
6 Transport and Storage.....	63
7 Delivery.....	64
8 Setting and Adjustment Possibilities.....	64
8.1 Back.....	65
8.2 Armrest.....	66
8.3 Control Panel.....	66
8.4 Footrest	68
9 Getting Into and Out of the Wheelchair.....	69
9.1 From the Side	69
9.2 From the Front.....	69
10 Putting into Operation	70
10.1 Switching On.....	71
10.2 "Drive" Function.....	72
10.3 Additional "Mode" Function.....	72
10.4 Lighting.....	72
10.5 Brake Release	73
10.6 Activation of Drive-away Lock	75
10.7 Horn	76
10.8 Battery Indicator	76
10.9 Charging	76
10.10 Batteries	79
10.11 Fuse	81
11 Status and Error Analysis	82
11.1 Warning.....	82
11.2 Errors	83
11.3 Defect / Failure	84
12 Options.....	84
12.1 Lift Seat.....	85
12.2 Electric Seat Tilt.....	87
12.3 Spring-mounted Front Wheels	87
12.4 Front Wheel Swivel Lock	88
12.5 Elevating Footrests	89
12.5.1 Height Adjustment	89
12.5.2 Elevating the Footrests.....	89
12.6 Electric Footrests	89
12.7 Electric Back Angle Adjustment of the Standard Seat.....	90
12.8 Recaro Seat	91

	Page
12.9 Control Panel.....	91
12.9.1 Lower Position Mounting.....	91
12.9.2 Swing-away Control Panel	92
12.10 Delta Control Panel	93
12.11 LCD	93
12.12 Blow-and-Suck Control.....	94
12.13 Push Button Module	96
12.14 Chin Control	96
12.15 External Mileage Indicator.....	97
12.16 Lap Belt.....	98
12.17 Other Options.....	98
13 Maintenance and Cleaning	99
13.1 Maintenance.....	99
13.2 Defective Lighting	100
13.3 Changing Tyres.....	102
13.4 Cleaning	103
14 Disposal	103
15 Technical Data of the B600.....	104
16 Technical Data of the Battery Charger	105
17 Terms of Extended Warranty	106
17.1 Subject of the Extended Warranty.....	106
17.2 Scope of the Extended Warranty	106
17.3 Conditions for Laying Claims	106
17.4 Exclusion of the Extended Warranty	107
17.5 Liability	107
17.6 Incidental Provisions	107

1 General Information

The present Instructions for Use describe the manipulation of the B600 power wheelchair. They provide the user and the attendants with the necessary knowledge regarding function, maintenance and disposal of the wheelchair. They contain all information which is needed for a safe manipulation as well as possible causes and help for elimination in case of malfunctions.

These Instructions for Use have been prepared on the basis of the

DIN EN 62079 standard "Preparation of Instruction Manuals - Construction, Contents and Representation".

The B600 power wheelchair has been constructed in accordance with the currently valid technical rules and is safe in operation. For the safety, harmonized standards, regulations for prevention of accidents and directives apply.

The B600 power wheelchair is a quality product offering versatile use in everyday life for the mobility of disabled people.

The following special features are the characteristics of the B600:

- Individual adaptations are possible through options and custom fabrication with modular components (chassis, seating system, control unit, accessories).
- Easy servicing due to good and clear accessibility of all component groups.

Before using the B600 power wheelchair, the user or attendant must read the Instructions for Use thoroughly, especially the section "Safety Instructions".

The design, as described in these Instructions for Use, is subject to technical alterations without notice.

2 Declaration of Conformity

Otto Bock HealthCare GmbH as manufacturer with sole responsibility declares that the B600 power wheelchair conforms to the requirements of the European Directive for Medical Products 93/42/EEC.

Development, construction and design of this product fully meet the requirements regarding technical safety laid down in the following standards:

- DIN EN 12182 Technical products for disabled
"General Requirements and Methods of Testing"
- DIN EN 12184 Power wheelchairs and mobiles and corresponding battery
chargers
"Requirements and Methods of Testing"

3 Safety Instructions

3.1 Warning Symbols

Note!
Information regarding manipulation of the product.

Attention!
Warning messages regarding possible technical damages.

Danger!
Warning messages regarding possible risks of accident or injury.

Note!
Information regarding environment protection.

3.2 General Safety Instructions

Danger!
To avoid potentially dangerous situations such as tipping, it is necessary that you become familiar with your new power wheelchair on level ground first.

Danger!
When using the wheelchair on inclined surfaces ensure the automatic brakes are operating correctly.

Danger!
Be especially careful with any naked flames and cigarettes, as the back upholstery and seat cushion could catch fire.

**Danger!**

When using your wheelchair in the public traffic, you are subjected to observe the traffic regulations according to the Road Traffic Act.

**Attention!**

When getting into or out of the power wheelchair, the wheelchair control must be switched off.

**Attention!**

Do not use the footrests as a stepstool when getting into or out of your wheelchair.

**Attention!**

Get to know with assistance from another person how the power wheelchair reacts when the center of gravity is shifted, i.e. when driving on slopes or inclines or clearing obstacles like steps and curbs.

**Attention!**

Use your power wheelchair properly.
The critical obstacle height is 5 cm. Therefore, it is necessary to avoid driving over steps or curbs higher than 5 cm.

**Attention!**

Do not drive against obstacles (including steps, curbs or doors) without braking first.

**Attention!**

When using lifting platforms, be sure that the power wheelchair is standing centrically on the lifting platform area and that none of the wheelchair's components, such as the anti-tipper, is protruding into the danger area.

**Attention!**

When using the wheelchair on lifting platforms, elevators or in the underground, the wheelchair control must be switched off and the brake must not be released.

**Note!**

The maximum load for the B600 power wheelchair is 120 kg (265 lbs.). Optionally the B600 power wheelchair can be retrofitted for a maximum load of 180 kg (397 lbs.).

**Note!**

The driving characteristics of the B600 power wheelchair can be influenced by electromagnetic fields (mobile phones or other radiating devices). All mobile devices must be switched off when driving.

**Note!**

It is also possible that the B600 power wheelchair itself generates electromagnetic fields that might cause interference in other devices. Therefore, the control is to be switched off whenever you don't need it.

**Note!**

Ensure the tyres have sufficient tread depth and check that the pneumatic tyres are inflated to the correct tyre pressure (printed on the sidewall of the tyre).

**Note!**

Ensure you are visible in the dark. If possible, wear light clothes or clothing with reflectors. Be sure the reflectors installed at the side and the rear of the power wheelchair are clearly visible.

**Note!**

The B600 power wheelchair may be used within a temperature range from -25 °C to +50 °C.

**Note!**

The B600 is not suited for driving on very slippery (such as icy) surfaces.

**Note!**

Avoid driving on very coarse-grained surfaces (gravel, pebbles, shingle or scree).

Note!
The B600 power wheelchair has been approved for driving on inclines and slopes of up to 17 %.

Note!
During battery charging, the wheelchair control must be switched off.

Note!
For safety reasons, we recommend attaching a lap belt, which is available as accessory.

Note!
Never use a water hose or high-pressure cleaning apparatus for cleaning the wheelchair.

Note!
When the B600 is not used over extended periods of time or for dispatching the B600, the fuse must be removed.

4 Field of Application

The B600 power wheelchair is designed solely for people who are unable to walk or who have a walking impediment. The wheelchair can be controlled either by the patient or by an attendant with the use of special controls.

The B600 is a power wheelchair that is designed for indoor and outdoor use. It is sufficiently compact and manoeuvrable for the use indoors and can overcome some obstacles outdoors (category B of the EN 12184 standard).

The versatility of this wheelchair is due to its modular design with adaptable options and adjustment possibilities. Typical use would be for patients with walking impediments / inability due to:

- Palsies
- Paralysis
- Loss of limbs
- Defective / deformed limbs
- Joint contractures / defects
- Other diseases

The B600 power wheelchair was specially designed for users who are able to independently move in a power wheelchair.

Fitting Considerations:

- Body height and weight (maximal load is 120 kg, optionally 180 kg)
- Physical and psychological limitations
- Age of the user
- Home conditions and environment

5 Service

Service and repairs on the B600 power wheelchair may only be carried out by authorised dealers. Should any problems arise, please contact your wheelchair supplier.

Your authorised dealer is:

Note!

The B600 power wheelchair must be checked and serviced at least ONCE a year by an authorised dealer!

6 Transport and Storage

When transporting the B600 power wheelchair in a vehicle, ensure the wheelchair control is switched off and that the automatic brake has not been disengaged. The wheelchair must then be fully secured by means of tension straps.

The frame of the B600 power wheelchair has four fixing eyelets for attachment of the straps.

Using your Otto Bock product as a seat in motor vehicles

We recommend that, wherever and whenever possible, wheelchair users transfer to the seats installed in the motor vehicle and use the corresponding vehicle restraint systems, because this is the only way to ensure optimum protection of the passengers in case of an accident.

It is, however, possible to use your Otto Bock **B600** as a seat in a motor vehicle, provided that the Otto Bock safety elements and appropriate restraint systems are used. For more information, please refer to our instructions for use:

"Using your wheelchair/ mobility base with seating shell or buggy as a seat in motor vehicles", article number: 646D158

Warning!

Please read the 646D158 Instructions for Use before using your Otto Bock product as a seat in motor vehicles!

To reduce the folding size for transport of the wheelchair, fold down the backrest and remove the armrests and footrests (Fig. 1).

Fig. 1: Folding size

For further details please refer to the section "Setting and Adjustment Possibilities". The temperature during transport and storage must be between -10 °C and +40 °C.

**Note!**

When the B600 is not used over extended periods of time or for dispatching the B600, the fuse must be removed!

When driving onto ramps the swivel lock, which ensures effortless straight line forward movement of the B600, is very helpful.

This swivel lock is optionally available and can also be easily retrofitted.

7 Delivery

Upon delivery by the specialist dealer the B600 is ready for use.

All settings have been made in accordance with the indications on the order form or are made directly on site by the dealer.

The B600 is adapted to the personal requirements of the individual.

Before using the wheelchair for the first time, all components of the standard equipment (Fig. 2) must be checked for completeness.

- | | |
|---|-------------------------|
| 1 | Backrest |
| 2 | Armrest |
| 3 | Joystick & control unit |
| 4 | Seat cushion |
| 5 | Footplate |
| 6 | Brake release |
| 7 | Motor & rear wheel |
| 8 | Anti-tippers |

Fig. 2: Standard equipment

The functions of the individual components must be checked according to the instructions in the sections "Setting and Adjustment Possibilities" and "Putting into Operation". Should any problems arise please refer to the section "Status and Error Analysis".

The tool kit provided contains one set of Allen wrenches, sizes 3 mm, 4 mm, 5 mm, 6 mm, and one open-end 13 mm spanner.

8 Setting and Adjustment Possibilities

The B600 power wheelchair offers numerous setting possibilities.

The seat height, seat width and seat angle have been set in accordance with the indications on the customer order form and may only be changed by the dealer.

Setting of the rear wheel spring system which is included in the standard equipment must be carried out by the dealer.

The telescoping front frame of the B600 power wheelchair allows the wheelbase to be adapted to the seat in use (Fig. 3).

Since this setting involves a shift of the center of gravity it may only be carried out by an authorised dealer.

However the user can vary the back angle, change the armrest height and lower leg length and disassemble the footrests and armrests themselves.

Fig. 3: Telescoping front frame

8.1 Back

The back can be adjusted to four different angle positions by pulling the release strap which is located at the base of the back.

Once the back is in the desired position, release the strap to engage the locks (Fig. 4).

Fig. 4: Release strap

8.2 Armrest

Remove the armrests by loosening the thumb screws and lifting the armrest. The thumb screws are located at the lower end of the armrest holders (Fig. 5).

Fig. 5: Removal of the armrests

The armrest height is changed by loosening the screw at the upper end of the armrest holders (Fig. 6).

Ensure the screw is snugly re-tightened after making any adjustments.

Fig. 6: Height adjustment of the armrests

8.3 Control Panel

For adjusting the position of the control panel on the armrest, loosen the attachment screws on the underside of the armrest.

Slide the control panel backwards or forwards to position it to the desired arm length (Fig. 7).

Fig. 7: Adjustment of the control panel

This function is very helpful when using channel forearm pads and hand pads, which are described in the Wheelchair Accessories Catalogue. If the control panel rail is too long, the protruding part can simply be cut off with a saw.

After loosening the connection plug, the control panel can also be removed completely from the armrest. To do this, turn the lock ring of the plug to the left up to its end position (Fig. 8).

Fig. 8: Loosening of the connection plug

Pull off the plug. When remounting the control panel, make sure to turn the lock ring to the right until it locks in place.

For further adjustment possibilities of the control panel please refer to the section "Options".

8.4 Footrest

To disassemble the footrest, the calf band must be unhooked. Release the footrest lock (Fig. 9), then swing away the footrest either inwards or outwards. When in this position, the footrest can be pulled upwards and removed.

Fig. 9: Footrest lock

Take care not to get your fingers jammed when swinging the footrest away or back in place.

For remounting, relocate the footrest in either the inner or outer positions and swing it back into position so the lock engages.

Reattach the calf band to the holding device.

To adjust the height of the footplates loosen the screw at the footrest hanger and adjust the height to suit your lower leg length and compensate for the seat cushion height used. The footplates adjustment must not be extended by more than 160 mm below the hanger (Fig. 10).

Fig. 10: Adjustment of the footrest hanger

**Attention!**

Ensure all screws and nuts are firmly re-tightened after making any adjustments.

9 Getting Into and Out of the Wheelchair

The B600 power wheelchair has been designed as a modular system.

The armrests and footrests can be removed easily and thus allow easy access for transfer from the side or from the front.

9.1 From the Side

Move the B600 power wheelchair as close as possible to the seat.

For this purpose, remove the armrests on the side where the user will get into the wheelchair (Fig. 11) and also remove the footrest if necessary.

Now the user can slide onto the seat.

Use of a transfer board will make transfer easier.

The lighting installation can be swung down for transfer.

Fig. 11: Removal of armrest

9.2 From the Front

After flipping up the two footplates, the user can get into and out of the wheelchair from the front (Fig. 12).

The space for getting into/out of the wheelchair can be increased by swinging the footrests to the side (Fig. 13).

The footrests can also be removed completely.

With the assistance of an attendant or by means of a transfer lifter, the user can get into or out of the wheelchair from the front. Use of a rotation plate is also possible.

Fig. 12: Flipping up the footplates

Fig. 13: Swinging away the footrests

10 Putting into Operation

To control the B600 power wheelchair, the Modular Control System is used. The Modular Control System comprises a control panel for entering the driving instructions by the user and for displaying the current state as well as a controller, which uses the input data to control the motors and other electric functions. Data transmission is realized via a bus system. This modular design allows to connect other modules and input devices such as chin control, suck-and-blow control or electric seat functions. The specialist dealer can program the control and thus adapt the wheelchair control to the individual user's requirements.

10.1 Switching On

Note!

Upon delivery, the main fuse of the B600 power wheelchair is not connected. It is contained in the protective covering on the control panel. Before using the B600 power wheelchair, this fuse must be inserted into the fuse holder at the rear end of the battery case (see section "Putting into Operation / Fuse").

Press button 1 to switch on the wheelchair control (Fig. 14).

The indicator 4 shows the remaining battery capacity.

The indicator 10 will show 1-4 lit LEDs depending on the selected speed level. When the battery indicator flashes or the LEDs continually light up one after another, this means a system error has occurred.

When the control unit is not used for some time, the wheelchair is switched off automatically. The control can be switched off at any time using button 1; if this occurs whilst driving the wheelchair will be braked immediately.

Fig. 14: Control panel

Item	Description
1	On/Off
2	Horn
3	Light
4	Battery indicator
5	Joystick
6	Direction indicators, right and left
7	Warning flasher
8	Indicator "Drive or additional function"
9	Increase speed level / Additional function
10	Speed level indicator

10.2 "Drive" Function

After switching on, the control is in the speed level which had been selected last. Button 9 is used to increase the speed level. After speed level 4, the control switches back to speed level 1. For driving use the joystick 5. The further the joystick is moved away from the mid-position, the faster the B600 power wheelchair will drive in this direction. The maximum speed with full deflection of the joystick is dependent on the selected speed level. The characteristics for speed, acceleration and deceleration can be adapted by the dealer to meet the individual user's requirements.

As soon as the joystick is released, the wheelchair is braked automatically until it stops moving. When standing still, the mechanical brakes are automatically active so the B600 cannot move.

If the drive function is not required, the wheelchair control must be switched off to prevent uncontrolled driving through accidental joystick manipulation.

10.3 Additional "Mode" Function

When the button M is pressed and held for approx. 2 seconds, the system will change from driving mode to the additional functions. In the speed level indicator, each LED stands for a certain additional function. To change between the different functions, shortly press the button M or move the joystick to the right. To switch back to the driving mode, again hold the button M depressed for approx. 2 seconds.

10.4 Lighting

Standard equipment of the B600 power wheelchair includes a lighting installation.

One halogen light and one front direction indicator are attached to the front frame on the right and left side.

For user transfer to and from the side of the wheelchair, they can be folded away (with standard seat) (Fig. 15).

To protect the direction indicators against damage they are provided with rubber bearings. Rear direction indicators and rear lights are integrated into the rear covering (Fig. 16).

The control panel is used to switch on and off the warning flasher, lights and direction indicators.

For changing defective lamps, please refer to section "Maintenance and Cleaning".

Fig. 15: Folding away the front lights

Fig. 16: Rear lights

10.5 Brake Release

In case of failure of the control unit or insufficient battery capacity, it is possible to push the B600 power wheelchair. To do this, the brake must be released using the mechanical release mechanism. The brake release mechanism is located either on the right or left side of the wheelchair between the frame and seat.

For safety reasons you must first pull up the release bolt (Fig. 17).

! Note:

When in this position, the brake must still be engaged! Should the brake already be released and the pushing function be active at this stage, please contact your wheelchair supplier.

! Note:

Repairs and adjustments on the brake may only be carried out by specialist staff.

! Danger

Wrong adjustments can lead to brake function failure and danger to life!

Fig. 17: Brake release bolt

Now the brake release lever has been released and can be pressed forwards and downwards fully (Fig. 18).

In this position, the control will recognize that the brake has been released and automatically deactivate the driving function.

Fig. 18: Brake release lever

Attention!

If the control shows no error signal when the brake is released, the brake setting must be checked by the specialist dealer.

Attention!

After turning the brake release lever into the pushing position, all brake systems are switched off.

To reactivate the brake, pull the lever upwards until the release bolt engages.
To reactivate the driving function switch off and then switch on again the wheelchair control.

10.6 Activation of Drive-away Lock

Activation via hand control (hand control with light, control panel B)

To activate the drive-away lock press the warning flasher button longer than 5 seconds. The control will acknowledge this with a short acoustic signal and then switch off automatically.

Activation via control stick

Since there is no warning flasher button in the system with this system configuration, use the following new menu item for activation.

Drive mode

- Drive 1
- Drive 2
- Drive 3
- Drive 4
- Drive-away lock

Additional functions

...

After selecting the entry "Drive-away lock" (USER) the control will acknowledge this with an acoustic signal and then switch off automatically.

Switching the System on with activated Drive-away Lock

System without LCD

After switching on and the LED test, the speed level LEDs will quickly light up one after another (LED 1, 2, 3, 4, 3, 2, 1, 2, 3, 4, 3, ... with only one LED active at a time). Now move the joystick to the front end position until a beep signal occurs after approx. 3 seconds. Next, move the joystick to the rear until again a beep signal occurs after approx. 3 seconds. Now the drive-away lock is deactivated and driving is possible. If the joystick is not moved correctly, the drive-away lock remains activated, even if the joystick is moved once again and this time correctly. After 10 seconds the control switches off again automatically.

System with LCD

After switching on, an information window appears on the LCD. The following example illustrates the kind of display. For the rest, the system is analog to the version without LCD (see above).


Fig. 19: Drive-away lock

10.7 Horn

Press button 2 to sound the horn.

10.8 Battery Indicator

The battery indicator displays different functions:

- Readiness for operation of the control unit
- Remaining battery capacity
- Status and error messages

When in driving mode, the indicator is lit on a steady basis and shows the remaining battery capacity. The indicator consists of 7 segments, each segment representing approximately 14 % of the total capacity. With a distance range of approx. 35 km on even ground at constant speed, each segment corresponds to a distance range of approx. 5 km.

Information regarding correct use of the battery indicator:

Battery is charged

Charge the battery, if possible

Battery charging is urgently required

Immediately after switching on, the indicator shows the battery capacity status saved before the wheelchair was switched off the last time.

The exact state of the remaining capacity will be displayed after driving for a short moment.

10.9 Charging

The batteries' capacity determines the distance range of the B600 power wheelchair.

Many factors have an influence on the batteries' capacity. Besides the temperature, the age of the batteries and the time and extent of driving, the way you charge the batteries has considerable effects on the capacity and thus on the distance range.

Driving over a longer period of time with only the red LEDs of the battery indicator lit will result in battery discharge and thus battery damage. In addition, there is the risk that the B600 may stop due to zero battery capacity and bring the user into a dangerous situation.

If the B600 power wheelchair is used daily, we recommend connecting the battery charger over night to always obtain full capacity and thus the maximal distance range.

If the B600 is not used for a longer period of time, a charging cycle must be carried out once a week to maintain the batteries' capacity.

Note!
The warranty is void in case of damages that result from a totally discharged battery!

Note!
When the B600 is not used over extended periods of time, disconnect the main fuse or charge the B600 power wheelchair once a week.

Note!
Only the battery charger provided by Otto Bock HealthCare GmbH may be used for charging.

Before connecting the battery charger ensure the control is switched off, now connect the battery charger to the charging receptacle (Fig. 20). The charging receptacle is integrated in the left or right lateral covering. Next, connect the battery charger to a wall socket.

The indications on the type plate of the battery charger must be identical to the voltage of the mains supply of your country.

Charging will start automatically.

The current charging state is indicated by the LEDs of the battery charger.

Battery Charger

The battery charger is designed for maintenance-free (gel) and low-maintenance (acid) lead batteries with a nominal voltage of 24 V and a nominal capacity of C5= 50..80 Ah. For optimal charging of both kinds of batteries two charging characteristic curves are saved in the charger. Upon delivery of the wheelchair, the battery charger's characteristic curve will be factory-set in accordance with the kind of battery the wheelchair is equipped with. If the battery charger is to be used with another wheelchair or if new batteries are installed in the wheelchair, the characteristic curve must be checked. The kind of battery can be changed by means of a switch located on the underside of the charger (covered by a green adhesive label) and is indicated by a flashing signal of the green LED at the beginning of the charging process.

Attention!
Wrong setting of the characteristic curve can result in permanent damage to the battery!


Fig. 20: Charging receptacle

**Note!**

During battery charging, the control must be switched off.
The complete charging current will be fed into the battery.

The battery charger features 3 LED indicator lamps:

Yellow LED is lit		Battery is charging
Green LED if flashing	1 time at the beginning of the charging process	Characteristic curve for low-maintenance lead batteries (acid)
	2 times at the beginning of the charging process	Characteristic curve for maintenance-free lead batteries (gel)
	Permanent flashing:	Battery is charged to approx. 80%
Green LED is lit		Battery is fully charged
Red LED is lit		Defective battery; battery is missing; charging time has been exceeded
No LED is lit		No mains supply

The battery charger features a programmed recharging phase.
After approx. 8 hours of full charging of a previously discharged battery, the battery charger can remain connected without having to fear overcharging or damage to the battery. When the charging process is complete, first disconnect the battery charger from the mains supply and then disconnect the battery charger from the charging receptacle.

**Note!**

The rubber cap must be re-inserted into the charging receptacle to protect it from humidity.

Switching the wheelchair control on again will make the B600 ready for driving.

**Note!**

During battery charging, explosive gases can develop. For this reason, provide for sufficient ventilation when using the battery charger in closed rooms.

Do not smoke when being near the battery charger or handling the batteries!

Avoid fire and sparking!

Do not cover the vent slots!

- Use the battery charger only within the indicated temperature range and humidity limits (see section "Technical Data of the Battery Charger").
- Always place the battery charger with its rubber feet on a level surface.
- Protect the battery charger from direct sunlight, since this would additionally heat up the charger.
- Avoid dust and dirt from affecting the function of the battery charger.
- To clean the charger, use a dry cloth.

10.10 Batteries

The standard version of the B600 power wheelchair includes two 12 V lead acid batteries with a capacity of 60 Ah.

During the charging process, the water acid mixture is converted into gas for a short time. The liquid level must be checked on a monthly basis and be topped up if required.

For maintenance purposes and for mounting/dismounting the batteries, the B600 power wheelchair features a battery drawer (Fig. 21).

Fig. 21: Battery drawer

Push the release bolt, which is located in the recess below the bumper (Fig. 22), upwards. Now you can easily pull out the drawer to the rear.

Fig. 22: Release bolt for battery drawer

To remove the cover of the battery case undo the snap lock of the safety strap (Fig. 23).

Fig. 23: Snap lock of the safety strap

Note!

Before making any maintenance work on the batteries, please read the warning instructions of the battery's manufacturer thoroughly!

Once the cover is removed, you have free access to the batteries from above (Fig. 24).


Fig. 24: Battery cover

Note!

When checking the battery's filling level of the B600 power wheelchair, make sure the B600 is standing on level ground!

Unscrew the cover caps of the individual chambers using a large screwdriver to see the battery's liquid level on the acid level markings. If the liquid is below this measure, top up with distilled water.

Maintenance-free lead gel batteries can be used as an option.

**Attention!**

When closing the battery drawer make sure that the release bolt correctly snaps into place. The release bolt must engage in the drive unit sustainer in the recess provided for this purpose! If correctly locked, the bolt is visible and clearly tangible on the underside of the drive unit sustainer.

Fig. 25: Release bolt engaged

10.11 Fuse

The 80 A fuse is located at the rear end of the battery case.

To replace the fuse, undo the lateral snap lock to open the cover of the fuse holder (Fig. 26).

After replacing the fuse, make sure that the cover properly snaps into place.

**Attention!**

During any maintenance work on the B600 where the battery cover is open, the fuse must always be removed for safety reasons.

Fig. 26: Fuse holder

11 Status and Error Analysis

The MCS control has different indication categories for the status and error analysis. Depending on the effects on the system, the control differentiates the following indications:

- Warning
- Error
- Defect
- Failure

In the following tables, the indication codes with the corresponding error groups as well as the possible causes and measures are described.


11.1 Warning

A warning indicates a status or malfunction of individual components.

Components without error are not restricted in their function.

This means, for example, that the light function remains intact in the case that one direction indicator has failed.

Display: **The LEDs flash 2 times per second for a duration of 4 seconds.**
In addition, 2 acoustic signals are emitted.
This process is repeated every minute or happens once when the system is started.



Indicator	Error group	Possible cause	Possible Measure
	Battery warning	Battery voltage too high or too low	Charge or discharge depending on the battery indicator
	Temperature warning	Overheated impulse control (controller) or motor	Reduce driving speed / let cool down
	Impulse control warning	Motor regulation problem	Restart
	Multifunction button or On/Off button	Connection error / Button defect	Check connection cable / buttons
	Driving is restricted	Certain seat function has been applied (lift seat)	Bring seat in normal position
	Input / output module warning	Motor is not connected, overloaded motor, output (e.g. defective direction indicator)	Check connection cables / lamp contacts
	Battery charger blocking	Battery charger is connected	Disconnect the battery charger after charging

11.2 Errors

An error affects the function of the MCS control. The system is not fully operational until the error has been eliminated.

Display: **Flashing signals of the LEDs: 4 x / second**

Acoustic signal: 2 x / second.

Indicator	Error group	Possible cause	Possible measure
	Battery error	Battery is discharged or overcharged	Charge or discharge depending on the battery indicator
	Right drive motor or brake	Right drive motor or brake are not connected, poor contact or cable break	Check motor cable and plug connection on the controller
	Left drive motor or brake	Left drive motor or brake are not connected, poor contact or cable break	Check motor cable and plug connection on the controller
	Error of input module	Joystick is not in mid-position, not calibrated during system start, defective input module	Check the mechanical mid-position and that the joystick moves easily / restart
	Error of output module	Defective connection to the output module / defective module	Check plug contact for firm seat / restart
	Error of impulse control	Error in the controller	Check plug contacts for firm seat / restart
	Communication error	Module is not connected, defective connection	Check plug contacts for firm seat / restart

11.3 Defect / Failure

A defect represents a serious failure of a system component.

After elimination of the error, the system will be re-activated by re-starting it.

A failure represents the most serious possible error status, which results in an immediate emergency stop of the system.

Display: **The LEDs continually light up one after another accompanied by an acoustic signal**

Indicator	Error group	Possible cause	Possible measure
	Incorrect system configuration	On/Off module is missing or defective, two equal modules are competing, wrong input device configured / connected	Connect the input devices, check the plug connection of the modules, restart
	Higher control device defective	Lost of the position of the joystick	Check the mechanical mid-position and if the joystick moves easily / restart

If you don't manage to directly eliminate the occurring errors with the described measures, the specialist dealer has the possibility to read the exact error code with the hand programming device and make a better system analysis.

The errors which occurred before are saved as error history and give information on all malfunctions which ever have occurred.

This will help to determine the further procedure of repair of the system.

12 Options

The B600 power wheelchair has been designed as a modular system.

The MCS control offers the possibility to use various electrical additional functions as well as specially adapted input devices.

In addition, other accessory components can be added to the wheelchair.

The complete range of accessory components is listed on the Order Form and in the Wheelchair Accessories Catalogue.

Some important options are described in the following.

12.1 Lift Seat

The B600 power wheelchair can be optionally equipped with a lift seat. This option allows to lift the seat by up to 40 cm by means of a motor drive. To control the seat lift function, either the joystick is used in the "additional function" mode or the push-button module is used. The driving function remains available even with lifted seat. Because of the lower stability of the power wheelchair against tipping, the maximum speed is reduced as soon as the lift seat leaves the lowest position. On the control panel this is indicated by flashing of the "Drive" LED. If an LCD is connected, a "snail symbol" appears to indicate the creep speed.

Fig. 27:
Snail symbol indicates creep speed

To provide safety when using the lift seat function, the following instructions must be observed:

Danger!
Use the lift seat function on level ground only.

Danger!
The driving operation with reduced speed is only to be used for manoeuvring.

Attention!
Maximum load capacity of the lift seat is 100 kg (220 lbs).

Danger!
Use the lift seat function only with vertically adjusted backrest.

Danger!
In the section between seat frame and wheelchair frame there are pinch points which are caused by the construction. This risk of getting pinched must be explained to attendants.

Attention!
Ensure that there are no objects in the lifting range.

**Danger!**

If the wheelchair does not switch to creep speed as soon as the lift seat function is activated, immediately contact your wheelchair supplier to eliminate the fault. Until fault repair, the power wheelchair may only be used with the seat in its lowest position.

**Attention!**

The use of the lift seat function (lifting and lowering) is restricted to a duration of 6 minutes per hour.

**Danger!**

When the mounting position of the lift seat device or of the patient seat is changed, there is the risk of tipping of the wheelchair.

**Danger!**

When retrofitting or changing the seat lift device, the instructions given in the Service Instructions must be observed.

12.2 Electric Seat Tilt

The electric seat tilt function allows seat tilt of up to 30° from the horizontal (Fig. 28). The integrated shift of the center of gravity improves the stability with maximal seat tilt. To control the seat tilt, the joystick is used in the "additional function" mode.

Fig. 28: Electric seat tilt

Attention!

Be particularly careful when driving on inclines with the seat in tilt position, because the stability against tipping is reduced!

12.3 Spring-mounted Front Wheels

Standard equipment of the B600 power wheelchair includes spring-mounted rear wheels. To increase driving comfort, the B600 offers the possibility to be equipped with spring-mounted front wheels, which provide for better adherence to the ground, especially when driving on uneven terrain (Fig. 29).

Fig. 29: Spring-mounted front wheels

12.4 Front Wheel Swivel Lock

In most cases when driving onto ramps or into elevators the wheelchair requires to be driven precisely in a straight line.

The front wheel swivel lock option of the B600 power wheelchair offers this possibility by locking the front wheels in the forward position.

This prevents the wheelchair from turning.

Slightly rotate the folding lever on the front frame to the side.

The folding lever will jump out of its normal mid-position (Fig. 30).

Fig. 30: Engaging the swivel lock

A bolt is released from the frame and engages in the front caster fork as soon as the straight forward position is reached. The wheelchair can no longer turn.

Now the B600 can drive straight forward or backward.

To release: Turn the folding lever back to the mid-position where it will disengage the bolt from the front caster fork (Fig. 31).

Fig. 31: Disengaging the swivel lock

12.5 Elevating Footrests

12.5.1 Height Adjustment

By loosening the screw at the footrest hanger, the footrest height can be adjusted to suit the individual lower leg length and compensate for the seat cushion height you use. There are five threaded holes in the footrest hanger for realization of the five positions of height adjustment.

Note!
Snugly re-tighten all screws and nuts.

12.5.2 Elevating the Footrests

Either apply a light pressure to the footplate and raise the leg or raise the footrest by hand to the required height.

After you have reached the desired height, the footrests will engage under weight bearing in one of the ratchet points.

To release the footrest lock, slightly lift the footrest and depress the lever located on the outside of upper footrest pivot point.

When you let go of the release lever, the footrest will engage in the next ratchet position.

12.6 Electric Footrests

To avoid permanent pressure loads, the B600 offers the possibility for mounting electric footrests (Fig. 32).

In the control's additional function menu, you can operate the right and left footrest separately as well as both at the same time.

Fig. 32: Electric footrests

12.7 Electric Back Angle Adjustment of the Standard Seat

It is possible to equip the standard seat with an electric back angle adjustment (Fig. 33).

**Attention!**

Be particularly careful when driving on inclines with the back inclined to the rear, because the stability against tipping is reduced!

Fig. 33: Back angle adjustment

To obtain a low folding height for transport, release the cross bolt at the lower end of the drive unit by pushing the lever upward (Fig. 34). Remove the side panels and fold down the back part onto the seat bottom. For unfolding, unfold the back and press the drive unit onto the lock to engage the bolt.

Fig. 34: Release bolt for back angle adjustment

**Note!**

Make sure that the backrest firmly locks in place!

12.8 Recaro Seat

The B600 power wheelchair can be equipped with different Recaro seat models.

The turning handle on the right or left side of the backrest is used to adjust the back angle.

To fold down the backrest, a release handle is used, which is located laterally at the backrest. By pulling the release handle upward, the backrest can be folded down to the front.

The Recaro seat can be equipped with an electric back angle adjustment. To control the electric back angle adjustment, the menu "Additional functions" of the wheelchair control is used.

The Recaro seat can easily be removed from the chassis.

Way of proceeding:

- Pull at the release strap which is located below the anterior part of the seat bottom.
- Loosen the two locking bolts which lock the seat to the right and left seat frame.
- To release the posterior seat attachment from the retainer bushings, slightly tilt the seat to the rear and slide back a little.
- Remove the Recaro seat.
- Now you have free access for any maintenance work.
- For remounting, place the posterior edge of the Recaro seat on the end of the seat frame.
- Slide the seat forwards until the posterior seat attachment hangs in the retainer bushings.
- Tilt the seat to the front until the locking bolts engage with the seat frame.

**Attention!**

Make sure when mounting the Recaro seat, that the seat is securely attached to the retainer bushings and that the anterior locking bolts have engaged up to the key ring!

12.9 Control Panel

12.9.1 Lower Position Mounting

A special control panel holder can be mounted as an option allowing the control panel to be mounted and adjusted downwards.

To do this, loosen the attachment screw (Fig. 35) and adapt the height of the control panel and joystick height to the armrest.

**Attention!**

Make sure to snugly re-tighten the screws and nuts after all adjustments!

Fig. 35: Control panel, lower position mounting

12.9.2 Swing-away Control Panel

To allow the user to drive the B600 power wheelchair closer to a table or even below the edge of a table, the control panel can be swung to the side with a special control panel holder.

Way of proceeding:

- By applying slight pressure, press the control panel holder to the side to release the pivot element.
- The control panel can now be swung away to the back up to the armrest (Fig. 37).
- When bringing the control panel back to its original position, the pivot element engages automatically.

Fig. 36: Release lever

Fig. 37: Control panel, swing-away to the side

12.10 Delta Control Panel

The two-piece Delta control panel can be mounted to allow angle adjustment of the control panel.

Loosen by approx. 1 turn the adjustment screw which is located below the black cover cap.

Now you can change the angle of the indicator / push button module.

Firmly re-tighten the screw after making adjustments.

The Delta control panel can be used together with the option "control for attendant and table-top control".

12.11 LCD

When special controls such as chin control or blow-and-suck control are used, the LCD module is used as indicator module.

When using the control panel, the LCD module is very helpful (Fig. 38).

In the driving mode the display shows the speed level, speed and the distance covered.

To increase the speed level, shortly press the mode button on the control panel.

The joystick is used for the drive function.

To change from driving mode to menu mode hold the mode button depressed for a while (Fig. 39). Here you can change the speed level and control additional functions such as lights and horn.

Fig. 38: LCD module

Fig. 39: Change from driving mode to menu mode

Menu control:

Joystick to the front	>>	go up in the menu
Joystick to the rear	>>	go down in the menu
Joystick to the right	>>	select in the menu
Joystick to the left	>>	one level upward in the menu
Mode button	>>	go back to driving mode

When different input devices are used, the wheelchair will be controlled by the module which was used to switch the B600 on.

12.12 Blow-and-Suck Control

When the blow-and-suck control is used, the drive and additional functions are controlled by different combinations of blow and suck instructions.

In addition, the multifunction switch (MFS) is used for moving within the menu (Fig. 40).

With slow driving speed, the steering position is displayed on the LCD module.

Fig. 40: Blow-and-suck control

There are two types of instructions:

- Short sucking or blowing
- Long sucking or blowing

Menu control

Control ON:	long pressure on MFS
Control OFF:	long pressure on MFS
Go up:	2 x short blowing
Go down:	2 x short sucking
To the right / select:	1 x short blowing
To the left:	1 x short sucking
One menu level upward:	2 x short sucking / short pressure on MFS

Continued mode

The wheelchair moves until a new instruction is given.

Forward travel:	1 x short and 1 x long blowing
Reverse motion:	1 x short and 1 x long sucking
Right turn during driving:	1 x short blowing
Left turn during driving:	1 x short sucking
Right turn directly after starting:	1 x long blowing
Left turn directly after starting:	1 x long sucking
Emergency stop:	long pressure on MFS
Change speed level (1,2...,5,1,2..):	2 x short blowing
Leave driving mode:	2 x short sucking

Un-continued mode

The wheelchair moves as long as the instruction is given.

Forward travel:	1 x short and 1 x long blowing
Reverse motion:	1 x short and 1 x long sucking
Right turn:	1 x long blowing
Left turn:	1 x long sucking
Emergency stop:	Push on MFS / long sucking or blowing
Leave driving mode:	2 x short sucking

12.13 Push Button Module

The push button module offers the possibility to control electric additional functions during normal driving. 5 pairs of push buttons allow to control the seat tilt, electric back adjustment, lift function, as well as the right and left footrest.

Fig. 41: Push button module

12.14 Chin Control

The chin control allows the user to control the wheelchair functions with a small joystick which can be positioned directly at the user's chin by means of a swivel arm.

A separate multifunction switch is used to control all required switch functions (Fig. 42). The multifunction switch has three possible positions.

The I/O button (on top) is used to switch the control on and off. The mode button (below) can be used (as with the normal control panel) to increase the speed level and – when depressed for a while – to leave the driving mode.

The swivel button (right) is used to swivel the joystick arm.

As long as this button is being pressed, the joystick arm moves up to the user.

When pressed again, the joystick arm will change direction and move away from the user.

Another swivel switch is attached to the rear side of the headrest to also allow the attendant to control this function.

The joystick's functions are similar to those of the standard joystick.

Fig. 42: Chin control

**Attention!**

When a special control is used, the specialist dealer must instruct the user in detail on how to use the control!

Other special controls such as LED-push button control, 5-button keyboard control, Buddy-Button control and environment control are available in conjunction with the LCD module.

12.15 External Mileage Indicator

As an optional extra you can have fitted to the protective bow of the control panel an external mileage indicator, which includes a speedometer, a day's mileage indicator, a total mileage indicator and a time indicator (Fig. 43).

Fig. 43: Mileage indicator

If the display shows m/h in the left upper corner, this means the speed indicator is active. Using the yellow button you can switch to the day's mileage indicator which shows three flashing arrows on the left side of the display.

To reset the day's mileage indicator to zero, depress and hold the button for 3 seconds. Shortly press the button again to display the total mileage indicator.

Press again to see the time display.

To set the time, depress and hold the button for 3 seconds.

12: and **24:** will appear in alternating manner.

To activate the 24-hour display, press the button whilst the **24:** is displayed, to activate the 12-hour display, press the button whilst the **12:** is displayed.

Now the time is displayed and the digits are automatically advanced one after the other.

To set the time, depress the button until the desired digits are shown on the display.

When the total mileage indicator is active you can change between m/h and km/h by pressing and holding the yellow button for 3 seconds. For the mileage measurement the wheel circumference needs to be entered here.

To calculate the wheel circumference use the following formula:

Wheel size 14"**Inches in cm** **14" x 25.4 = 355.6 mm****U = d x π** **35.56 mm x 3.1416 = value to be entered**

12.16 Lap Belt

The B600 power wheelchair can be equipped with a lap belt (Fig. 44).

The lap belt serves exclusively as a safety belt when driving the power wheelchair.

! Important: The lap belt must never be used as a safety belt for the user when being transported in the wheelchair in a motor vehicle.

Fig. 44: Lap belt

12.17 Other Options

- **Head support fixture** (Fig. 45)

For the attachment of headrests (see Accessories Catalogue)

Fig. 45: Head support fixture

-
- **Curb Climbing Assist**
For clearing curbs and steps having a maximal height of 10 cm
 - **Puncture Resistant Tyres**
Solid rubber tyres
 - **Anatomically Shaped Seating Systems**
Different models of "Recaro" seats are available for greater comfort and more adaptation possibilities
 - **Belt Systems**
Lap belt, chest belt, and four-way chest strap kit from Schroth
 - **Seating Shell Interface**
For the adaption of special seating shells; control panel holder for seating shells
 - **Armrest Accessories**
Special adapters for the armrests (see Accessories Catalogue)
 - **Attendant Control**
Separate control panel at the backrest
 - **Joystick Accessories**
Fork for tetraplegics, STICK S80, soft ball, golf ball, flexible control stick shaft
 - **Protective Bow for Control Panel**
Metal bow for protecting the control panel against shocks and damage
 - **Crutch Holder**
 - **Folding Rearview Mirror**
 - **Tray**
Attachable tray top
 - **Accessories Case, Accessories Bag, Pocket for mobile phone**

Other optional add-on components are available on the Order Form and in the Accessories Catalogue.

13 Maintenance and Cleaning

13.1 Maintenance

The function of the B600 power wheelchair must be checked before each use. The items listed in the following table must be checked by the user at the indicated intervals. Failure to carry out these simple checks may lead to problems arising that could invalidate the warranty.

	To be checked	daily	weekly	monthly
Tyres	Tyre air pressure (printed on the side wall of the tyre)			X
	Sufficient tread depth			X
Batteries	Liquid / acid level (not with gel batteries)			X
Lighting	Damage of the housing Electric function	X	X	
Electronics	Control without error message	>> each time before driving <<		
	Battery charger doesn't show any error messages on the LEDs		X	
Brake	Manipulate brake lever with the control switched on	X		
Brake	Adjustment of the brake release lever			X

**Attention!**

When the brake lever is manipulated, the wheelchair control emits an error signal on the control panel. If this function is not ensured, the specialist dealer must make a corresponding setting!

If you have any problems carrying out this maintenance work, please contact an authorised dealer (section "Service").

13.2 Defective Lighting

If a housing or lamp is defective, the corresponding components can be ordered from your specialist dealer.

To change the halogen lamp of the front light, slightly press the black engaging lever at the lower end of the lamp. Now the pane can be folded to the front and removed (Fig. 46).

Take the lamp by the lock ring, pull it out and replace it with a new one.

When mounting, make sure that the notch on the lamp socket is correctly positioned.

Fig. 46: Replacing the front light

For changing the rear direction indicators and the rear lights, unscrew the pane. Lightly press in the lamp and turn it to the left by approx. 90°, then remove it (Fig. 47). Insert the matching lamp and turn it to the right to provide function. Further information on the lamps used is contained in the section "Technical Data".

Note!

When re-inserting the panes make sure they are seated exactly on the housing and check that the screws are firmly tightened. Avoid that moisture can enter the lamps!

Fig. 47: Replacing rear direction indicators and rear lamps

13.3 Changing Tyres

Way of proceeding:

- For dismounting the rear wheels, loosen the Phillips head screw in the middle of the wheel and remove the wheel cover (Fig. 48).
- Loosen the screw of the wheel (Fig. 49) and pull off the wheel.
- Disassemble the two-piece rim by loosening the four Allen head screws.
- The defective tube is freely accessible and can be replaced.
- For dismounting the front tyres, loosen the screw of the axle and remove the axle.
- Again, disassemble the two-piece rim and replace the defective tube.

Fig. 48: Wheel cover

Fig. 49: Screw of the wheel

13.4 Cleaning

When cleaning the B600 power wheelchair, the electrical components must not come in contact with water.

Note!

- For cleaning the control panel, battery charger, armrest and seat upholstery use only a damp piece of cloth and a mild cleansing solution.
- For cleaning the backrest and seat cushion use a dry brush.
- For cleaning wheels and frame use a damp plastic brush.
- Please prevent the electronics, motors and batteries from direct contact with water.
- Do not use any scouring agents.
- Never use a water hose or high-pressure cleaning apparatus for cleaning the wheelchair.
- Please ensure your B600 power wheelchair is checked and serviced at least ONCE a year by an authorised dealer!

14 Disposal

For disposal of the B600 power wheelchair, the disposal regulations of the respective country shall apply.

15 Technical Data of the B600

Measurements and weights

Seat width:		38 - 42 cm or 43 - 48 cm
Seat depth:		38 - 46 cm or 42 - 50 cm
Seat height:		45 - 60 cm
Armrest height:		24 - 36 cm
Armrest length:		26 cm
Lower leg length:		25 - 34 cm or 35 - 44 cm or 45 - 54 cm
Back height:		45 or 55 cm
Back angle:		-9/1/11/21° or 0/10/20/30°
Overall width:		64.5 cm
Overall height:		103 cm
Overall length:		108 cm
Turning radius:		155 cm
Tyre size:	front	8" - 10"
	rear	14"
Air pressure:	front	printed on the sidewall of the tyre
Weight when empty:		115 kg
Maximum load capacity:		120 kg (optionally 180 kg)

Electrical installation

Operating voltage:		24 V
Batteries:	wet cell batteries	2 x 12 V, 60 Ah (5 h)
	gel batteries	2 x 12 V, 50 Ah (5 h)
Electronics:	model	MCS
		24 V DC max. 100 A
Lighting:	front direction indicator	H21W 12 V BAY9s
	front light	HMP 08 2.4 W, 6 V, PX13.5s
	rear direction indicator	C21W 12 V, BA15s
	rear light	C5W 6 V, BA15s
Battery charger:	model	MEG 2408
	input	230 V, 60/50 Hz
	output	24 V, DC/8 A
Fuse:		80 A

Driving data

Speed:	6 km/h / 10 km/h
Climbing ability:	17 %
Obstacle height that can be cleared:	5 cm
Distance range:	approx. 35 km

Technical modifications reserved.

16 Technical Data of the Battery Charger

Automatic battery charger with computer-controlled characteristic curve for 24 V lead acid batteries.

Power requirements:	230 V -10 %, +6 %
Mains frequency:	50 Hz ± 4 %
Protection class:	1 (protective conductor)
Charging connection:	24 VDC
Nominal charging current:	8 A
Residual ripple:	<1 %
Charging characteristic curve:	IUU with trickle charging phase, computer-controlled
Primary fuse:	T4G fuse link, not accessible from outside
Secondary fuse:	Electronic, reversible reverse battery protection, short circuit-proof, idling-proof, overheating protection
Protective system:	IP 21
Ambient temperature:	-10 °C to + 40 °C
Indicator:	2 LEDs
Weight:	1.7 kg
Dimensions (WxHxD):	105 x 65 x 205 mm

17 Terms of Extended Warranty

17.1 Subject of the Extended Warranty

- This extended warranty applies to wheelchairs.
- The extended warranty comprises claims of the medical supply store / service company against Otto Bock HealthCare GmbH and does not affect the statutory warranty claims of the ultimate consumer towards the medical supply store or any other service company being responsible for the fitting of the patient.

17.2 Scope of the Extended Warranty

- Otto Bock HealthCare guarantees, subject to the adherence to the conditions described in paragraph 3 and subject to the observance of the exclusions according to paragraph 4, that beyond the statutory warranty period, up to 4 years after the first use no construction faults, fabrication faults or material faults will occur on frame parts and cross braces of manual wheelchairs fabricated in series production and on frame parts of power wheelchairs fabricated in series production that considerably impair the wheelchair's usability.
- If 2 years after the first use construction faults, fabrication faults or material faults occur, Otto Bock HealthCare has the choice either to replace the entire wheelchair or to replace or repair defective parts of the frame and of the cross braces. This extended warranty does not entitle the medical supply store / service company to further claims for reduction, cancellation of the contract or damage.

17.3 Conditions for Laying Claims

The medical supply store / service company can lay claims under the extended warranty, only if

- it refers to first use,
- the wheelchair has been employed by a healthcare provider authorised by the health-insurance system,
- the healthcare provider authorised by the health-insurance system passes the complete fitting on along with a report of defects,
- for the construction of custom fabrications made from the products fabricated in series production to adapt them to the individual conditions of a patient exclusively medical products with CE-mark are used and these products are used in conformity with the functional purpose appropriate for such medical products,

- modifications on a wheelchair do not go beyond the adaptation intended by Otto Bock HealthCare (e.g. adjusting the legrest length) and the use of the accessories recommended by Otto Bock HealthCare (modular system),
- the one who repairs (adapts) medical products uses exclusively spare parts and add-on parts authorised by the manufacturer (original spare parts) and works according to the provisions given by the manufacturer (see instructions for use).

17.4 Exclusion of the Extended Warranty

The extended warranty is not effective, if Otto Bock HealthCare proves that

- the conditions for laying claims according to paragraph 3 are not complied with,
- the diminution of the wheelchair's usability is due to unprofessional employment of the wheelchair especially retrofits which do not conform to the manufacturer's description,
- the diminution of the wheelchair's usability is due to normal wear and tear especially of motors, battery chargers and batteries whose usability is usually limited to one year,
- the diminution of the wheelchair's usability is due to incorrect storage, transport or inappropriate as well as unprofessional use and storage,
- the diminution of the wheelchair's usability is due to a change of the patient's physical constitution, e.g. considerable increase in weight,
- the diminution of the wheelchair's usability is due to force majeure.

17.5 Liability

- For the undertaking of performance under this extended warranty, Otto Bock HealthCare is liable only for slight negligence in case unessential obligations have been violated; Otto Bock HealthCare assumes no liability for intentional acts of persons performing extended warranty service on our behalf. In addition, damages claims are restricted to the damage foreseeable at the beginning of the undertaking of performance.

17.6 Incidental Provisions

- Otto Bock HealthCare is allowed to destroy replaced parts three weeks after return unless the patient or his/her paying party (health insurance company) expressly object.
- Place of fulfillment of performance under this extended warranty is Duderstadt.
- This extended warranty is governed by German civil and commercial law; place of jurisdiction is Duderstadt.

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Otto Bock HealthCare GmbH

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Otto Bock HealthCare GmbH has been certified by the German Society for the Certification of Quality Assurance Systems (DQS)
in accordance with DIN EN ISO 9001 standard, reg. no. 779 (management system)